



FLOUT® DOSING CHAMBER- INTERNAL VENT



FLOUT® DOSING CHAMBER- EXTERNAL VENT



FLOUT® Dosing Chamber RECOMMENDED MINIMUM DIMENSIONS for 3 inch /75mm OUTLET

Flout® body is 8 inches/200mm wide Minimum width of chamber 12"/300mm

Multiple outlet Flout® body is 8 inches/200mm per section. Outlet pipes 8 inches/200mm ON CENTER. Total minimum width of chamber is number of outlets x 8 inches/200mm plus 4 inches/100mm.

Avoid placing INLET over the Flout®. Provide a baffle to 4 inches/100mm off floor if necessary.

Flout® and Floating Outlet® are registered trademarks for a dosing apparatus manufactured under US Patent No. 5290434 by :

Rissy Plastics, LLC 350 Cedar Lane Torrington, CT 06790 860-482-1645 rissyplastics.com